Agilent 2100 Bioanalyzer System

Specifications

Sample throughput	11/ 12 Samples per chip in 30-60mins (Assay dependent)
Sizing Range	50 to 7000 bp
Sample volume	1 μΙ
Application	 RNA quality check – the industrial standard for RNA analysis, offering total RNA and mRNA data with RNA Integrity Number (RIN), and small RNA analysis. DNA size and quantity – high resolution separation and quantification of DNA Protein analysis – determination of the quality and purity of proteins from Coomassie down to silver stain sensitivity.

© BioCRF